

**FOR IMMEDIATE RELEASE**

Contact: Gary Deger / Steve Wesson

# Press Release

## New Thermal Desorption System can handle a Variety of Sorbent Tube Sizes

CDS Analytical has re-imagined the Thermal Desorption Autosampler. Announcing the Dynatherm TDA 9300, the first Cryogen-free thermal desorption autosampler on the market able to handle any size thermal desorption tube. The autosampler is perfect for replacing old Tekmar 6016 units and Nutech VOST systems, automating chemical agent monitoring or implementing product emissions testing in an R & D lab.

Our new 18-tube autosampler is able to handle tubes from 3" to 71/2" long, with ODs of 3mm to 16mm, complete with computer control, automatic leak checking, rapid heat rates, and no need for bulky cooling systems.



SEE US AT ARAB LAB STAND NO. 208

465 Limestone Road, P.O. Box 277, Oxford, PA 19363-0277 USA  
610-932-3636 Fax: 610-932-4158 e-mail: [gdeger@cdsanalytical.com](mailto:gdeger@cdsanalytical.com) Web: [www.cdsanalytical.com](http://www.cdsanalytical.com)

**CDS** Analytical, Inc.