



### **pfm Slide 4003 E Sliding microtome (electronic)**

The newly designed electronic sliding microtome pfm Slide 4003 E is a high - performance instrument with adjustable retraction for all paraffin sectioning applications in routine use, research and industry.

The object movement as well as the trimming function and retraction are controlled via a noise - free stepper motor.

All functions can be operated manually or automatically and can be activated, changed or deactivated. The display indicates the current settings.

The fine- and coarse-adjustment of the object is operated by a button that is integrated into the front lateral working area of the microtome. This setting controls the advance of the object to the knife edge.

As well as a memory function and a section counter fixed reference points in the lower and upper section area facilitate setting-up the samples.



## **pfm Waterbath 1000**

The modern new developed pfm Waterbath 1000 offers you the possibility of stretching paraffin sections to be used in histological routine, research and industry laboratories.

### **Features**

- Housing: light anodised aluminium
- Internal tank : black anodised aluminium
- Easy to clean
- Fast and precise adjustment of the water temperature
- Save heating system to avoid overheating and damage even in case of empty waterbath
- Water protected power supply
- Water level control
- Blue marking to show where to empty the water
- Illumination of the water (white or blue colour optional)

### **For further information please contact:**

pfm Produkte für die Medizin AG  
 Wankelstrasse 60  
 50996 Köln  
 Germany  
 Tel. 02236/9641-0, Fax -20  
 E-Mail: [erwin.schneider@pfm-ag.de](mailto:erwin.schneider@pfm-ag.de)  
[www.pfm-ag.de](http://www.pfm-ag.de)